DENOMINAZIONE / PRODUCT NAME / NOM DU PRODUIT / PRODUKTBEZEICHNUNG / DENOMINACIÓN / NOME DO PRODUTO

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Eurosteril® sterilization pouches and rolls are used for the sterilization of medical-surgical devices and of materials for medical use which are especially suitable for steam or EtO GAS. They are made by 2 layers: one made of high resistance medical paper and the other made of blue polyester/polypropylene. They are available as prefolded and self-adhesive rolls and pouches of various sizes in order to satisfy different needs of packaging. Sealing temperature between 170 °C and 200 °C. The seal is made by two parallel strings which are separated by a groove for rolls and sterilization pouches with a width of 5 cm. The seal is made by three parallel strings which are separated by two grooves for rolls and sterilization pouches with other dimensions. On the paper layer two chemical indicators of process are printed: one for steam and the other for ethylene-oxyde, both indicating that the sterilization process has been carried out.

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PROPERTIES	UNITS	STANDARDS OR METHODS	TYPICAL (TOLERANCE)
THICKNESS	µm	-	42 (± 10%)
BASE WEIGHT	g/m ²	-	45 (± 10%)
BREAKING LOAD	kg/mm ²	ASTM D 882	25 MD (± 5) 25 TD (± 5)
ELONGATION TO BREAK	%	ASTM D 882	130 MD (± 40) 130 TD (± 40)
WATER VAPOUR PERMEABILITY	g/m ² /24h	ASTM F1249 38°C- 90% RH	≥ 8 (± 10%)
OXYGEN PERMEABILITY	cc/m ² /24h	ASTM D 3985 - 23°C - 0% RH	≥ 127 (± 10%)

**CARATTERISTICHE
PRESTAZIONALI /
PERFORMANCE
CHARACTERISTICS /
CARACTÉRISTIQUES DE
PERFORMANCE /
LEISTUNGSINHALTE /
CARACTERÍSTICAS DE
RENDIMIENTO /
CARACTERÍSTICAS DE
DESEMPENHO**

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Eurosteril® pouches and rolls are intended to contain medical devices, resist the sterilization process allowing the sterilization of the contents, allow the preservation of the contents in sterile conditions.

Eurosteril® pouches and reels are compatible with the following sterilization treatments:

- STEAM
- EtO (Ethylene Oxide)

Pouches and reels are equipped with Type 1 process indicators (EN ISO 11140-1) for Ethylene Oxide and Vapor, printed with water-based ink. These provide an indication for the user as they identify whether the barrier system, together with its contents, have already undergone one of the sterilization treatments indicated above. Indicators follow the criterion:

- Steam sterilization = STEAM indicator from pink to brown
- Sterilization with ethylene oxide = EO indicator from blue to yellow / brown

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The sealing on the 4th side for pouches, and on the 3rd and 4th side for reels, can be carried out with all thermosealing machines suitable for the construction materials of Eurosteril® products. It is advisable to refer to the International Standard EN 868-5 for the correct setting parameters (temperature, pressure, speed) in order to guarantee the desired minimum resistance according to the sterilization method employed.

The melting temperature of polypropylene is 175°C.

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Do not re-use. Do not use if package is damaged.

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Storage: sterilized pouches and rolls must be handled as less as possible, kept possibly in closed cabinets in a cleaned environment, at a temperature and humidity under control. If maintained as described above, at a temperature between 8 and 35 °C and a relative humidity of between 30 and 70%, Eurosteril packaging systems will keep the contents sterile for a maximum of 6 months and, in any case, no later than the indicated expired date (table 1 DIN 58953-8:2010).

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Medical devices must be disposed of in accordance with medical procedures adopted by health care institutions and organizations and local, state or federal laws and regulations.

Sterile Barrier System devices of the Eurosteril® series are comparable to urban waste and therefore they must be disposed of as urban waste (mixed, non-recyclable).

The PET plastic film has been processed in such a way as to obtain a stratification that makes it non-recyclable. It must therefore be disposed of as non-recyclable waste.

If the operator decouples film from paper, the paper can be recycled in separate collection.

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