




# Product Data Sheet 2368

Rev 1  
Effective date: 01.05.09  
Page 1 of 1  
Authorisation: CO1502

**Product Description** A single use test pack for Euronda E9 inspection® 18/24 Type B pre-vacuum STEAM sterilizers working at 134°C for 3.5 minutes

It consists of a number of specially selected sheets of barrier material and a chemical indicator sheet. The whole pack is wrapped in steam sterilizable paper. The indicator sheet inside the pack is sensitive to Time, Steam and Temperature (TST).

A shipper case contains 10 Bowie Dick tests and includes instruction leaflet, lot and expiration data.



**Specification** Colour changes from yellow to completely blue / purple if rapid and even steam penetration has occurred at 134°C for a holding period of up to 3.5 minutes.

**Applicable standards** EN 867-5:2001 Class B

**Product Size** 124mm x 124mm

<b>Packaging</b>	Primary case	10 Test Packs	Size	260mmx170mmx145mm	Weight:	3kg Gross
	Contents:					

**Recommended Storage** Cool, dark, dry conditions before and after use: 0 - 30°C 30 – 60 % RH

**Shelf Life** 36 months from date of manufacture

**Active Components** Each indicator, which is the indicating medium, contains 0.7mg dyes and 0.7mg reagents. None of the substances used in the formulation are known carcinogens, nor do they contain any heavy metals.

The component chemicals are potentially irritant. Being fully encapsulated, this should present no occasion for contact with the chemical components. If however, contact is made due to breakdown of the encapsulation, wash the skin with soap and water.

If ingested wash out the mouth thoroughly and give plenty of water to drink.

None of the components are of animal origin, hence to the best of our knowledge is free of TSE (Transmissible Spongiform Encephalopathy) agents.

**Safe Usage** The indicator and packaging contain no added heavy metals, no known carcinogens, and no added rubber latex. The indicator is considered safe when used under normal conditions.

The TST Single Use Bowie Dick Test is designed to be used in high vacuum steam autoclaves. When used in such devices no inherent hazards are likely to present themselves.

When used as intended , the indicator does not release any known to be toxic in sufficient quantities to cause a health hazard, during or after the sterilization process for which it is designated, in accordance with ISO 11140-1:2005, 5.9

**Disposal** Treat unwanted indicators and packaging in the same way as normal paper waste.

**Albert Browne Ltd, Chancery House,  
190 Waterside Road, Hamilton Industrial Park,  
Leicester, LE5 1 QZ, UK  
A Subsidiary of STERIS Corporation  
Tel: +44(0)116 2768636**

