Euronda | Alle®

EN INTENDED USE
Devices intended to convey the irrigation solution to the drill and/or to the surgical area to be treated during dental and maxillofacial surgery.
USER
The device is intended for use by qualified medical, technical and specialist trained personnel (ie. surgeon, dentist).
PATIENT TARGET GROUP
The intended target population is patients regardless of the gender, including pregnant and lactating women. Not for use on children below 2 years of age or below 10 kg of weight.

WARNINGS
• The use of the device is contraindicated in patients with sensitivity or allergies to PVC, ABS, Silicone.
• The content shall be presented following the appropriate aseptic presentation. Responsibility for aseptic presentation lies within professional personnel.
• Do not reuse the device. The reuse could lead to cross contamination.
• Always ensure correct operating conditions (ie. absence of leaks) and that sufficient and adequate coolant is delivered. Otherwise replace the device.
• After use, protect the spike with the cap provided.
• Report any serious incident that has occurred to the manufacturer and to the competent authority.

ASEPTIC PRESENTATION PROCEDURE
Prerequisites
• Personnel: 2 interacting people are required:
 a. a scrubbed person, wearing sterile gloves, face mask and gown as per the local applicable procedure;
 b. a second person to open the sterile barrier and to present the device. Clothing must comply with the provisions of the local applicable procedure.
• Area: adequate space, sufficient room, with an appropriately controlled environment and sufficient light to grant clear vision.
• Check for the integrity of the sterile barrier. Content is sterile unless the sterile barrier has been opened or damaged. If it is damaged or unintentionally opened before use, replace the device.
• Check the use-by-date to make sure that the device has not expired.

Presentation
• The person opening the sterile barrier shall be positioned close to the scrubbed personnel.
• Avoid involuntary contacts.
• The person opening the sterile barrier shall grab the 2 flaps on the side of the sterile barrier

where the symbol  is positioned and slowly peel the 2 flaps apart following the direction specified by the arrow. Refer to the visual description.
The device shall be presented to the scrubbed person making sure that there is no contact with any unsterile areas during the transfer.
The scrubbed person shall place the device in the sterile field.

SPECIFICATIONS OF USE
• Always consult the use of the instructions for use provided with the dental surgical motor related to the connection with the irrigation set for machinery.
• Connection procedure shall be conducted as follow:
1. Remove the cap from the spike.
2. Insert the spike into the bag. The spike has an air filter with an open-close-system, which allows for regulation of the flow rate.
3. Where necessary, use the roller to adjust the flow rate.
4. Insert the elastomeric section into the peristaltic pump, paying close attention to the flow direction.
5. If equipped with a tap, connect the end extension and adjust the direction of the flow.
6. Connect the final part of the irrigation set for machinery to the handpiece. The irrigation set for machinery offers the following options:
A. One-way internal irrigation
B. One-way external irrigation
C. Two-way internal and external irrigation, using the Y-connector provided.
7. Secure the irrigation set for machinery to the handpiece tube using the clips provided.

DISPOSAL
When the device is no longer needed, protect the spike with the cap provided. Once used, the device is contaminated by biological fluids. Follow your local unit and country-specific laws, directives, standards and guidelines for disposal.
DESCRIPTION OF THE SYMBOLS

	BLISTER PEELING DIRECTION Indicates the peeling direction to correctly open the blister.
	NUMBER OF PIECES PER BOX Indicates the number of units in each consumer pack.

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