**Composition and Description**  
Each 100 ml contains:  
Sodium Chloride: 0,9 g  
Water for injection q.s.: 100 ml  
Electrolytes:  
Sodium: 154 mEq/L  
Chloride: 154 mEq/L  
Irrigation solution is a sterile and non-pyrogenic solution for irrigation in PP bags of 250 mL, 500 mL or 1000 mL volumes.

**Indications**  
- Only for irrigation in medical and surgical applications.  
- Indicated for irrigation, cleaning of body surfaces and cavities, including eyes.  
- Also used for dilution of concentrated irrigation and cleaning solutions.

**Contraindications**  
- Not for use in surgical procedures related to the heart or central circulatory system.  
- Not for use in surgical procedures that directly contact the central nervous system.  
- Not for use intravenously.  
- Not for use in applications where electrocautery is used.

**Warnings**  
- Consult "instructions for use" before usage.  
- Do not use for intravenous injection.  
- Do not use for electrocauter applications.  
- Do not use if the solution is not clear, contains particles or if the bag is damaged.  
- Microwave heating of the solution is not permitted.  
- Aseptic technique is essential during application.  
- Product is for single use only. Discard unused portions of solution.  
- After the container is opened, the solution should be used promptly to minimize the possibility of bacterial growth or pyrogen formation.  
- Used products should be treated as medical waste.  
- Solution should not be used after the expiry date given on the label.  
- Solution should be used with care in patients with congestive heart failure or severe kidney impairment.

**User and Population**  
- There is no age or gender discrimination for intended use. Product shall be used by healthcare professionals only.

**Side Effects**  
- There can be adverse reactions if patient is sensitive to raw materials.  
- Excessive volume during irrigation of closed cavities like bladder may cause undue distention of tissues.  
- Lack of aseptic technique may cause site infection which may spread to other sites.  
- Absorption of large amounts of liquid can cause adverse reactions.  
- These adverse reactions caused by sodium or chloride amounts of solution are not expected during routine treatment. These adverse reactions frequency is not known.

**Dosage and Usage:**  
- Dosage to be applied will result from the type of administration. It can be dispensed by pouring into the area to clean, or through irrigation.

**Storage**  
- Store at room temperature between 15°C - 30°C  
- Keep this product out of the reach and sight of children.  
- Do not use this product after the expiry date.  
- Do not dispose of expired or remained products in the environment! Follow the local regulations for medical waste management.

**Manufacturer**  
**POLIFARMA I.L.A.C.SAN.VE.T.I.C.A.S.**  
Vakfir OSB Mah, Sanayi Cad. No:22/1, 59930 Ergene / Tekirdağ / TÜRKİYE  
www.polifarma.com.tr  
**Distributor**  
**EURODRA S.P.A.**  
Via dell'Artigianato 7, 36030 Montebelluna Peralto, Vicenza / ITALY  
www.eurodras.com

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