

## STERILE TRANSPARENT DRAPE COVER WITH ADHESIVE

Sterile drape cover disposable for medical devices manufactured of transparent polyethylene film.

- Low density polyethylene (LDPE) transparent film waterproof
- Adhesive tape doesn't leave any residue and sticks well over the skin

**TECHNICAL FEATURES:** Dimensional tolerance  $\pm 5\%$

### PRODUCTION:

The product is produced and packaged in a special **clean room with a controlled environment**

### CLEAN ROOM ENVIRONMENT:

ISO class 8 - UNI-EN ISO 14644-1 standard.

Constant temperature and humidity environment maintained by ventilation and filtering systems that also allows to keep the bioburden at low levels.

### PACKAGING:

Sterile drape, packed in bag for sterilization with peel-open opening and sterilization indication, packed in cardboard boxes.

### STERILIZATION:

Sterilization is carried out with ethylene oxide; the sterilization process is regularly validated and kept under constant control by routine analysis; and the residue of the sterilizing gas is kept within the limits established by the current legislation. Sterilization is performed by a specialized laboratory.

The sterilisation process is validated in accordance with the harmonised European standard EN 550 and the ISO 11135 standard, as well as with F.U.I. and F.E.

### SHELF LIFE :

The validity of the sterile device is 59 months from the sterilization date (MONTH-YEAR). The expiration date indicated refers to the device in its intact and correctly stored package.

**QUALITY CONTROL:** ALLE has implemented a Quality Management System which complies with the following standards:

- **EN ISO 13485**

### GENERAL PRECAUTIONS AND WARNINGS:

The device must be used only by qualified and trained personnel. Read the information reported on the label.

### STORAGE:

- The items are stored in the original cardboard boxes, in a closed and clean environment protected from sunlight
- The integral sales package is guaranteed impermeable to direct and indirect light sources.

### PRODUCT DISPOSAL:

The device listed in this document can be disposed of as non-hazardous waste in authorized plants, including incineration, in accordance with current regulations.

However, if the device used during surgical procedures in the operating room is exposed to contact with the patient's biological fluids, it must be treated and disposed of as hospital waste.

### LATEX AND PHTHALATE-FREE:

- The device does not contain natural rubber latex (Latex free)
- The device does not contain phthalates (PVC free)
- The device does not contain elements (adhesives, chemicals, etc.) that could cause allergies or any skin reaction.

**CND: T020199 (national classification of medical devices)**

**PRODUCT CLASSIFICATION: CLASS IS Medical Device**

## STERILE TRANSPARENT DRAPE COVER WITH ADHESIVE

DIMENSION DRAPE / ADHESIVE BAND	PACKAGING	CODE
cm 40x50 / cm 4x40	20 pcs	270227



TECHNICAL SPECIFICATIONS OF THE MATERIALS		
Low density polyethylene (LDPE) transparent film waterproof.		
MEDICAL-GRADE ADHESIVE		
	MEAN VALUES	TEST
THICKNESS	Adhesive 102 $\mu$ + Support 107 $\mu$ = Total 209 $\mu$	[TEST: ATM-T07/010]
COLOUR	Transparent adhesive/White support	
ADHESION	N/m 350 - 458 (90°C polyethylene adhesion)	[TEST: ATM-T04/097]
DETACHMENT FROM SUPPORT	g/m 1260	[ATM-T13/001]

