

TECHNICAL SHEET

OF PRODUCT

PRODUCT NAME	MIXING BOWLS
CODE	D1,D2,D3,D4,D14
MATERIAL	PVC PT65 NON-TOXIC
THERMAL PROPERTIES	A 30° TENDS TO SOFTEN
MECHANICAL, PHYSICAL PROPERTIES OR CHARACTERISTICS	SELF CLEANING LATEX FREE
CHEMICAL PROPERTIES	
SAFETY DATA SHEET	N/A
INSTRUCTION FOR USE	
CAUTIONS	NONE
RISKS ANALYSIS REF.D.L.46 24/2/97	N/N
CLASSIFICATION 93/42/ EEC	NO

N/A DOES NOT NECESSARY

N/N NOT NECESSARY

tel. +39 0445 329811

fax +39 0445 865246

www.euronda.cominfo@euronda.comEURONDA SPA
Via dell'Artigianato, 7
36030 Montebelluna PrecalcinoVicenza - Italy
Cap. Soc. € 175.440,00 i.v.Cod. Fisc. e P.IVA 00595740242
I.S.O. IT 00595740242Reg.Impr. di VI n. 00595740242
R.E.A. di Vicenza n. 139557