

<b>DISPOSABLE OVERSHOES IN CPE</b>	
<p>Disposable CPE overshoes in light blue color provides effective protection both in wet and dry conditions. Overshoes are used in the operating room as protection for the operator from infected liquids and to avoid pollution of the operating field or sterile and sanitized environments.</p> <ul style="list-style-type: none"> <li>• One-size overshoes, suitable for every type of foot</li> <li>• Overshoes with ankle elastic</li> <li>• Overshoes made of blue CPE (Chlorinated polyethylene)</li> <li>• Dimensions: length 36 cm, height 15 cm</li> </ul>	
<b>TECHNICAL FEATURES:</b> (Dimensional tolerance $\pm 10\%$ )	
<b>PRODUCTION:</b> The product is produced and packaged in a special <b>clean room with a controlled environment</b>	
<b>CLEAN ROOM WITH CONTROLLED ENVIRONMENT:</b> Class ISO 8 – Standard UNI-EN ISO 14644-1. Environment with constant temperature and humidity; by means of special ventilation and filtering systems, the presence of the bacterial load is kept controlled and at low levels.	
<b>PACKAGING:</b> Plastic Bag	
<b>STERILIZATION:</b> Articolo non sterile	
<b>QUALITY CONTROL:</b> The quality system adopted by ALLE complies with the following regulations: <ul style="list-style-type: none"> <li>• <b>EN ISO 13485</b></li> </ul> Checks are carried out at all stages of processing: entry, production, packaging.	
<b>GENERAL PRECAUTIONS AND WARNINGS</b> This product should only be used by qualified and trained personnel. Always read the label.	
<b>GENERAL PRECAUTIONS AND WARNINGS</b> This product should only be used by qualified and trained personnel. Always read the label. Storage is in the original cardboard boxes, indoors and in a clean environment away from sunlight	
<b>STORAGE:</b> <ul style="list-style-type: none"> <li>• The items are stored in the original cardboard boxes, in a closed and clean environment protected from sunlight</li> <li>• The integral sales package is guaranteed impermeable to direct and indirect light sources.</li> </ul>	
<b>PRODUCT DISPOSAL:</b> The device listed in this document can be disposed of as non-hazardous waste in authorized plants, including incineration, in accordance with current regulations. However, if the device used during surgical procedures in the operating room is exposed to contact with the patient's biological fluids, it must be treated and disposed of as hospital waste.	
<b>LATEX AND PHTHALATE-FREE:</b> <ul style="list-style-type: none"> <li>• The device does not contain natural rubber latex (Latex free)</li> <li>• The device does not contain phthalates (PVC free)</li> <li>• The device does not contain elements (adhesives, chemicals, etc.) that could cause allergies or any skin reaction.</li> </ul>	
<b>CND:</b> T0208 - OVERSHOES (National classification of medical devices IT)	
<b>PRODUCT CLASSIFICATION : CLASS I Medical device</b>	
<b>GMDN:</b> Not applicable	
<b>DISPOSABLE OVERSHOES IN CPE</b>	
<b>PACK</b>	<b>CODICE</b>
<b>100 PCS</b>	<b>270420</b>

