

Sandrigo 23/07/2026

Below the data relating to the correlation of concentration (ppm) of use and conductivity ($\mu\text{S}/\text{cm}$) in the washing cycle using Euronda Sanifizer Enzymatic and Euronda Sanifizer Neutralizer medical devices.

Device of test :

Termodisinfectore EUROSAFE60 (Ref S.N. ETM240012)

Cycle applied: cycle standard

	Concentration (%)	Electrical Differential conductivity ($\mu\text{S}/\text{cm}$)
Euronda Sanifizer Enzymatic	0.1	≤ 160
Euronda Sanifizer Neutralizer	0.1	≤ 318

* Data obtained by HPLC analysis

Euronda S.p.A.
Technical and R&D office