

PRODUCT DESCRIPTION	<p>A single use test pack for Euronda Type B pre-vacuum STEAM sterilizers working at 134°C for 3.5 minutes.</p> <p>It consists of a number of specially selected sheets of barrier material and a chemical indicator sheet. The whole pack is wrapped in steam sterilizable paper. The indicator sheet inside the pack is sensitive to Time, Steam and Temperature (TST).</p> <p>A shipper case contains 10 Bowie Dick tests and includes instruction leaflet, lot and expiration data.</p>
SPECIFICATION	<p>Colour changes from yellow to completely blue / purple if rapid and even steam penetration has occurred at 134°C for a holding period of up to 3.5 minutes.</p>
APPLICABLE STANDARDS	<p>EN 867-5:2001 Type B</p>
PRODUCT SIZE	<p>124mm x 124mm</p>
PACKAGING	<p>Primary Case Contents: 10 Test Packs</p> <p>Size: 260mm x 170mm x 145mm</p> <p>Weight: 3kg Gross</p>
RECOMMENDED STORAGE	<p>Cool, dark, dry conditions before and after use: 0 - 30°C 30 – 60 % RH</p>
SHELF LIFE	<p>36 months from date of manufacture</p>
ACTIVE COMPONENTS	<p>Each indicator contains 0.7mg dyes and 0.7mg reagents.</p>
SAFE USAGE	<p>The indicator is considered safe when used under normal conditions.</p> <p>The TST Single Use Bowie Dick Test is designed to be used in high vacuum steam autoclaves. When used in such devices no inherent hazards are likely to present themselves.</p>
DISPOSAL	<p>Treat unwanted indicators and packaging in the same way as normal paper waste.</p>

All technical information herein disclosed is provided to us by our raw material suppliers, under their own responsibility.