



SPECIFICATIONS OF USE

Consult the section of the Instruction for Use provided with the dental unit or surgical aspirator related to the dental unit or surgical aspirator to be used by the handpiece of the dental unit or surgical aspirator.

• Visually inspect the integrity of the product.

• Check the function of the device.

• Activate the suction system. The correct functioning implies the suction of water from the top of the canula without aspiration points, i.e., in the connection with the dental unit or surgical aspirator, using an adapter.

• After use, handle the device as specified in the REPROCESSING INSTRUCTIONS.

• The use of the device is considered safe if patient aspiration is not present.

• The mechanical action of the canula on the mouth must not cause irritation and superficial injuries.

• Devices are sterile unless delivered, before use, in a sealed container, which must be opened by the patient. Therefore, they must undergo complete REPROCESSING.

• Failure to comply with the REPROCESSING cycle may result in damage to the device.

• Before every use, the product must be visually checked for damage.

• In case of damage, the device must not be used.

• Clean and dry the device immediately after use.

• Dried or fixed residues can make cleaning more difficult.

• Report any serious incident that has occurred to the manufacturer.

SHIELD OF THE DEVICE

• Notified of all the specific protection parameters.

• If not notified after its last use, the device must be discarded as waste, since cross infection risks become higher.

• Alternatively, the sterilized device is to be disposed of according to the instructions for use.

• The use and routine monitoring of the device is the responsibility of the dental professional and personnel, achieves the desired results.

• The validation and routine monitoring of the device is the responsibility of the dental professional and personnel, achieves the desired results.

• The user is responsible for the correct classification of the medical products, the definition of the reprocessing procedure and the reprocessing times.

• **Reprocessing procedure:** standards, steps and specifications for the cleaning, disinfection and sterilization of medical products, as well as the application of the reprocessing times.

• **Post-treatment:** cleaning, disinfection and sterilization of the medical device ends with the documented clearance for storage and reuse.

• **Pre-treatment:** cleaning, disinfection and sterilization of the medical device before its first use.

• **Storage:** cleaning, disinfection and sterilization of the medical device before its first use and after every subsequent use.

• **Wear protective gloves:** wear protective gloves.

• **Use protective clothing:** wear protective clothing.

• **Use a face mask:** use a face mask.

• **Use protective brushes:** use protective brushes.

• **Never use brushes with metal bristles:** never use brushes with metal bristles.

• **Select a cleaning brush with nylon bristles to clean the external surfaces.**

• Select a round cleaning brush with nylon bristles for the cleaning of the inner surface.

• **Water temperature:** the water temperature (at least 20 °C) must be appropriate for the product.

• **Disinfectant:** a manual cleaning, intermediate rinsing and disinfection of the dental unit or surgical aspirator during a dental procedure, thus ensuring a dental surgery.

• **User:** the dental professional and the dental staff.

• **PATIENT/TARGET GROUP:** the dental professional and the dental staff.

• The aspiration tip is used in the context of an outpatient treatment, as well as in the context of an inpatient treatment, regardless of gender, including pregnant and non-pregnant women.

• **CHARACTERISTICS OF THE DEVICE**

EM19 L114mm; D16mm

EM21 L129mm; D16mm

EM40 ASPIRATOR TIP L200mm; D16mm

D

#

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

D

